

## Contents

<i>EX VIVO</i> RECEPTOR OCCUPANCY ASSAY .....	2
Animal dosing and tissue collection.....	2
Sectioning .....	2
Incubation in radiotracer .....	2
Phosphorimaging.....	2
Data analysis .....	2
<i>IN VIVO</i> RECEPTOR OCCUPANCY ASSAY .....	2
Animal dosing and tissue collection.....	2
Data analysis .....	3

## EX VIVO RECEPTOR OCCUPANCY ASSAY

### Animal dosing and tissue collection

Animals are administered the test article (oral, i.p., s.c or i.v) or vehicle alone. After allowing for tissue uptake and distribution, the animals are sacrificed and the brain or other organs removed. A blood sample is collected at the time of sacrifice for subsequent quantification of drug levels in the plasma.

### Sectioning

Frozen brains (or organ) from the drug-treated animals are trimmed with a razor blade and mounted in a cryostat chuck. Tissue sections are cut at a thickness of 20  $\mu\text{m}$  using a cryostat (Lieca® 3050S) and thaw mounted onto Superfrost® slides. Three consecutive sections are placed on each slide, with a total of two slides (6 sections) per brain.

### Incubation in radiotracer

Slides are placed horizontally in a humidified box and radioligand in assay buffer (1 ml) layered over each slide. Sections are incubated in the radioligand solution for 15 min at room temperature. The radioligand solution is then rapidly aspirated off and the slides immediately placed in ice-cold wash solution (three washes, 5 min each). Following the final wash the slides are dipped briefly in distilled water and dried under a stream of warm air.

### Phosphorimaging

Sections are placed over a multipurpose ( $^{125}\text{I}$ ) or tritium-sensitive ( $^3\text{H}$ ) phosphor screen together with autoradiographic standards. The screen is exposed for 1 - 5 days and then scanned on a Fujifilm® FLA-7000 scanner.

### Data analysis

Regions-of-interest (ROIs) are drawn over each section and radioactivity levels measured in units of psl/mm<sup>2</sup>. The psl/mm<sup>2</sup> values are converted to D.P.M/mm<sup>2</sup> by reference to the autoradiographic standards. A value for specific binding is generated by the subtraction of mean nonspecific binding from mean total binding for each brain or organ. Percent inhibition of specific binding is plotted against the drug dose or plasma or tissue drug concentration to determine receptor occupancy.

## IN VIVO RECEPTOR OCCUPANCY ASSAY

### Animal dosing and tissue collection

Animals are administered the test article (oral, i.p., s.c or i.v) or vehicle alone. After allowing for tissue uptake and distribution, the animals are placed in a restraint device and administered a radiolabeled tracer for the target receptor (in 0.2 ml saline) via a tail vein and returned to the home cage. 30 minutes later the animals are sacrificed and the brain or other organs removed.

A blood sample is collected at the time of sacrifice for subsequent quantification of drug levels in the plasma. Individual brain or tissue regions are dissected out and placed into preweighed vials and reweighed. Once all the tissue regions have been harvested the tissues are digested overnight in Solvable® (1 ml). Following digestion, scintillation cocktail (10 ml) is added to each vial and the vials counted in a liquid scintillation counter.

#### Data analysis

D.P.M values for each tissue region are converted to D.P.M/mg tissue and thence to % injected dose per gram. In vivo receptor binding of the radiotracer is quantified by determining the tissue radioactivity (D.P.M./mg) in a receptor-rich tissue area relative to that in a receptor-deficient, reference region. The latter comprises non-specific binding whereas the former comprises both receptor binding and non-specific binding. Receptor occupancy for the drug-treated animals is determined as a percentage of that in the vehicle-treated animals and the percentage receptor occupancy values plotted against plasma or tissue drug levels.